


Vertebral heart score— now with results trending and reference ranges

Gain even more insight with the latest enhancements to the vertebral heart score (VHS) tool:

- Trending graph—graph the last five VHS scores to monitor heart size over time.
 - Icons indicate manual versus automatic scoring.
 - Results and graph update in real time.
 - Click a data point to see the score and open the associated radiograph.
- Reference ranges—compare the patient's score with a reference range for breed and/or species.
- Color coding—**red** = out of range, **green** = in range, **yellow** = borderline.



The IDEXX Web PACS* vertebral heart score tool delivers valuable information for screening, tracking, and monitoring.

Look for new features and enhancements in every IDEXX Web PACS update.